ACCESSORIES


FOR MEDICAL PACKAGING



















Web









MANUFACTURER COMPANY SPECIALIZED IN MEDICAL SEALERS SINCE 1990.

We Are Making Firsts in the World

We are not only making things easier for you, we also pay attention to our carbon footprint in order to extend the life of our world and provide a better quality of life.

"Sleep Mode (Energy Saving System)" technology was applied by Gündem A.Ş. for the first time in our world in our "Sealers" and was exhibited at the Medica 2009 fair in Germany.

- Year 2006- "Device safety": Thanks to the Emergency Stop / Emergency Reverse / Emergency Forward technology, jammed pouches can be easily removed from the device.
- Year 2007- In order to prevent the produced pouches from falling apart and being collected easily; "Drawer System" for Workstations and Mobile Table.
- Year 2008- Sealing speed of 12m/min thanks to "Fast-Sealing" technology.
- Year 2009- Energy saving thanks to "Sleep Mode" technology -First in the World-
- Year 2009- Ability to print 2 data on 5cm pouches thanks to "Write-Matic" technology.
- Year 2010- Self-testing and reporting on compliance with ISO 11607-2 before each use thanks to "Validation-Test" technology.
- Year 2011- Sealing speed of 12m/min thanks to "Fast-Sealing" technology.
- Year 2012- "MultiLine Sealing" 14mm multi-barrier sealing
- Year 2013- Thanks to "Customized Teflon-Coating" technology, pouches that get stuck in the device and melt are prevented from sticking to the heaters.
- Year 2014- "Stand System" to free up space in the work area thanks to "Space-Saving" technology
- Year 2016- Production speed of 6000 pouches per hour with our "HPC High Production Capacity" 40GM Impulse Pouch Maker model
- Year 2018- Sealing speed of 16m/min thanks to "Fast-Sealing" technology
- Year 2019- "5 in 1" Rotary Pouch Maker device production with Internal Printer -First in the World-
- Year 2019- "Soft-Start" Soft start for heavy and large pouches
- Year 2020- Thanks to "Auto-Slider" technology, "Mode Changing" without removing the rolls from the device or changing the settings -First in the World-
- Year 2021- Ability to quickly drop to Tyvek pouch temperature thanks to "Auto-Cooling" technology -First in the World-
- Year 2022- Automatic adjustment of the gluing temperature of the device when pouch/roll type is selected, thanks to "Auto-Temp" technology
- Year 2023- Thanks to "PrintTo-Cut Side" technology, text can be printed on the cut sides of pouches on our Rotary Pouch Maker devices. -First in the World-
- Year 2023- Thanks to "Multi-Roll-Processing" technology, 8 rolls of 5cm or rolls of different widths up to 43cm can be processed simultaneously and in exactly the same length. -First in the World-
- Year 2023- Thanks to the "Roll is Over" technology, the device stops even if 1 roll runs out.
- Year 2023- Preventing unnecessary operation of the printer when there is no printer cartridge in the device, thanks to the "No-Ribbon" feature. -First in the World-and...
- Year We continue to work on R&D to offer you more innovations.

REGISTERED TRADEMARKS OF GÜNDEM A.Ş.


GÜNDEM® | ABGMat® | ABGPack® | ABGTape® | ABGIndicators®

CONTENT

Description	Short & Full Code	Page
Workstation XL	A105 GM1215WST	4
Workstation L+	A100 GM8817WST	4
Workstation L	A095 GM8815WST	5
Book Support for Workstations	A115 GM115BHL	5
Upper Shelf drawer for Workstations	A110 GM110CKB	6
Illuminated Magnifier for Workstations	A125 GM125MGL	6
Document Holder for Workstations	A120 GM120DHL	6
Ink Ribbon for Dot Matrix Printer	A175 GM175IRB & GM175IRBE	7
QR Code Scanner	A150 GM150QRR	7
Barcode Scanner	A135 GM135BRR	7
External Thermal Printer	A130 GM130ETP	8
Shelf Adhesive label roll for External Thermal Printer w/o Indicator	A131 GM131PLN	8
Shelf Adhesive label roll for External Thermal Printer with Steam & EO Indicators	A132 GM132PLM	8
Shelf Adhesive label roll for External Thermal Printer with H2O2 (Plasma) indicator	A133 GM133PLW	9
Resin ink ribbon for External Thermal Printer	A134 GM134PLR	9
Shelf Adhesive Thermal label rolls for Internal Thermal Printer	A160 GM160THL	9
Thermal paper rolls for Internal Thermal Printer	A155 GM155THP	10
Mobile Table Wide + Drawer System	A090 GM1000TFC	10
Mobile Table Wide	A085 GM1000TF	11
Mobile Table Narrow	A080 GM500TF	11
Autoclavable Basket Standard - 1 STU	A072 GM500KC	12

Description	Short & Full Code	Page
Autoclavable Basket Mini	A071 GM250KC	12
Manual Roll Cutter & Holder 50cm	A070 GM632CM	12
Wall Hanger and Rack-Cup Set for A070	A075 GM632CMA	13
Manual Roll Cutter 75cm	A065 GM632CC	13
Manual Roll Cutter 50cm	A060 GM631CC	13
Stand Wide 75cm	A055 GM632CCS	14
Stand Narrow 50cm	A050 GM631CCS	14
Roller Conveyor XL 111cm	A045 GM140RTXL	15
Roller Conveyor L 73cm	A040 GM140RTL	15
Roller Conveyor S 59cm	A035 GM140RT	15
Side Rack Set for Rotary Sealers & Manual Roll Cutters and Stands	A030 GM14WTA& A077 GM631CCA	16
Sorting Table for Automatic Impulse Pouch Makers	A025 GM450TC	16
Front Table XL for Rotary sealers with printer	A020 GM111TC	17
Front Table L for 06GM+ sealers	A018 GM040TC	17
Front Table S for narrow sealers	A010 GD030TC	17
Reel holder stand Large	A016 GD040RC	18
Reel holder stand Small	A015 GM030RC	18
Medical Packaging- Sterilization Reels	A999 Medipack	18
Sterilization Monitoring Indicators	A997 GM997STM	19
Identification Tapes for Surgical Instruments	A996 GM996IDT	19
Seal-Check-Test -Tyvek®	A170 GM170SCT	20
Seal-Check-Test - Standard	A165 GM165SCS	20








Technical Data					Su	ıitab	ole N	1ode	els				
This mobile workstation is built for Automatic Impulse Pouch Makers for packaging products, inspecting surgical instruments, etc. The Workstation includes a mobile base element, a large desk with upper and lower shelves, autoclavable 1STU baskets, Barcode reader and thermal printer racks.	W509	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W560	W580	+W590	05GD	04GD+
Contents 1 Dec Upper Storage Shelf Grp. (222-200-012-5.44)													
 1 Pcs. Upper Storage Shelf Grp. (223-300-012-SA4), 1 Pcs. Side Shelf Left Grp. Mini (223-300-015-SA5.2), 1 Pcs. Side Shelf Right Grp. (223-300-015-SA5), 1 Pcs. Study Table Grp. (223-300-012-SA2), 1 Pcs. Drawer System (With telescopic rail) (223-300-015-SA1.2.1), 1 Pcs. A072 Autoclavable Basket, 1 Pcs. Lower Storage Rack (223-300-012-SA1), 1 Pcs. Allen Key. Dimen.: 1245 x 821 x 1497mm (LxWxH), (Excl. device, basket and other accessories) Optional Accessories A115 Book Support, A110 Upper Shelf drawer, A125 Illuminated Magnifier, A120 Document Holder, A072 Autoclavable Basket (2 extra pieces can be installed). 													
Technical Data					Su	ıitab	ole N	1ode	els				
It is produced to support books, folders, etc. to be placed on the Top Shelf in all Workstations. It is not fixed to the top shelf.	A105	A100	A095										
• Entirely AISI 304 stainless steel. • Dimension: 290 x 170 x 200 mm (LxWxH), • Weight: 1.94 Kg													

External Thermal Printer A130 GM130ETP	Technical Data					Su	itak	ole M	1ode	els				
	This printer, which can work integrated with the device (depending on the features of the device), issues an adhesive label identifying the package content for each package or at the end of the specified target number of packages. These labels can be pasted on the packages or on the boxes in which the packages are placed. Labels are not erased in the autoclave. Ask for labels with indicators.	909	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W560	08GM	+W590	05GD	04GD+
Strong Resolution To Annual Mark Park	 Contents Printer, data and power cable, 1 roll label, 1 ink ribbon. Size: Size: 270 x 220 x 119 mm (LxWxH), Weight: 3.00 Kg 													
Shelf Adhesive Label Roll - for External Thermal Printer without Indicator A131 GM131PLN	Technical Data					Su	itak	ole M	1ode	els				
## Company	This Indicatorless label is produced for External Thermal Printer with code A130 GM130ETP. These labels can be sticked on the packages or on the boxes in which the packages are placed. Labels are not erased in the autoclave.													
Por Someone State of the South	• 1 roll label without indicator. • 500 labels per roll • Each label 100 x 100 mm													
Shelf Adhesive Label Roll - for External Thermal Printer with Steam & EO Indicators A132 GM132PLM	Technical Data					Su	itak	ole M	1ode	els				
POS. 50999000000000000000000000000000000000	This Steam and Ethylene Oxide indicator label is produced for External Thermal Printer with code A130 GM130ETP. These labels can be sticked on the packages or on the boxes in which the packages are placed. The indicators on the labels change color when processed in the sterilizer.	A130												
DO D	• 1 roll label with STEAM & EO process indicators according to the EN ISO 11140-1 • 500 labels per roll • Each label 100 x 100 mm													

Shelf Adhesive Label Roll - for External Thermal Printer with H2O2 (Plasma) indicator A133 GM133PLW	Technical Data					Su	ıitak	ole N	1ode	els				
POS. PEOSPESSANA POS. SOPSONO SORS WART 11112222233 (00) F 1222233 000000001 2	This Vaporized H2O2/Plasma indicator label is produced for External Thermal Printer with code A130 GM130ETP. These labels can be pasted on the packages or on the boxes in which the packages are placed. The indicator on the labels changes color when processed in the sterilizer. It does not contain cellulose and does not adversely affect sterilization.	A130												
POS PLASMA WASHINGTON 120 CONTRACT CONTRACT POS STANDARD POS STANDA	 Contents 1 roll label with H2O2/Plasma indicator 500 labels per roll Each label 100 x 100 mm 													
Resin Ink Ribbon - for External Thermal Printer A134 GM134PLR	Technical Data					Su	iitak	ole N	1ode	els				
1 Inch Core 80MM x 300M Resin Ribbon FOR Polyester Labels (Waterproof)	Printer ribbon. It is produced for External Thermal Printer with code A130 GM130ETP.	A130												
FOR Printing On Labels Of Up To 100mm Width FOR Printing On Labels Of Up To 100mm Width FOR Thermal Transfer Printer Which Uses 1" Core Ribbon	• 110 mm x 300 meters per roll													
Shelf Adhesive Thermal Label Rolls - for Internal Thermal Printer A160 GM160THL	Technical Data					Su	uitak	ole N	1ode	els				
POR SOPROBLEM SOR	This thermal label roll is made for the Internal Thermal Printer integrated into our devices. These labels can be sticked to bundled pouches and boxes to prevent disintegration of pouches produced on Pouch Maker devices and to facilitate stocktaking. It also prints out the Validation Reports of the devices. These labels can be sticked into reports. Thermal labels are not suitable for autoclaves.		20GM	46GM	45GM	40GM	30GM	20GM	10GM	W960	W580	+W590	05GD	04GD+
POS SOSSOCIOS SO	• 1 Box= 10 rolls • 7 meters per roll													

Thermal Paper Rolls - for Internal Thermal Printer A155 GM155THP	Technical Data					Su	uitak	ole N	1ode	els				
Gindeman ONAMONG OS:05	This paper roll is made for the Internal Thermal Printer integrated into our devices. These thermal papers can be attached to bundled pouches and boxes to prevent the pouches produced in Pouch Maker devices from falling apart and to facilitate stock counting. It also outputs the Validation Reports of the devices. Thermal papers are not suitable for autoclaves.	W509	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W960	W580	+W590	05GD	04GD+
Characters	• 1 Box= 10 rolls • 10 meters per roll													
Mobile Table Wide + Drawer System A090 GM1000TFC	Technical Data					Su	uitak	ole N	1ode	els				
	This mobile table is produced for packaging sterilized products, examining surgical instruments, etc. The mobile desk includes a mobile base element, a height-adjustable desk, Metal basket and Drawer System (Telescopic rail) with code "223-300-015-SA1.2.1" are optional.		20GM	46GM	45GM	40GM	30GM	20GM	10GM	MD60	08GM	+W590	05GD	10000
	 Contents Mobile table with four wheels, two of which have brakes, which is large enough to accommodate the device and the front table, and whose height can be easily adjusted by the user. This table is stainless steel oven painted and its top is AISI 304 stainless steel. It has a telescopic rail system for attaching an autoclave basket under the table. When not in use, this basket can be pushed under the table as a drawer. The pouches thrown by the Impulse Pouch Maker can be collected in the basket. A025 Grouping Table, which enables grouping/regular accumulation of pouches, is suitable for this table. It has a hook system on the sides to attach two 1 STU autoclave baskets when desired (optional). Dimensions suitable for Rotary Sealers with Printers and Impulse Pouch Makers and also for A055 Stand Size (Length x Width x Height): 603 X 780 X 694-994mm Weight: 20 Kg. Check out our narrow mobile table A080, suitable for other sealers. 													



Autoclavable Basket Standard- 1 STU A072 GM500KC	Technical Data					Sı	uitak	ole N	1ode	ls	
	This basket has 1 STU dimensions defined in EN ISO Standards. It can be easily attached and removed from its slots on Workstations and Mobile Tables and is resistant to autoclaves and sterilizers.	A105	A100	A095	A090	A085	A080				
	 Contents It has the size of 1 STU as defined in EN ISO Standards. All made of AISI 304 stainless steel. Dimension (L x W x H): 570 x 280 x 250 mm Weight: 2.45 Kg. 										
Autoclavable Basket Mini A071 GM250KC	Technical Data					Su	uitak	ole N	1ode	ls	
	This basket can be used on Workstations and Mobile Desks and is autoclave and sterilizer resistant.	A105	A100	A095							
	• All made of AISI 304 stainless steel. • 250X150X100mm • Weight: 0.40KG.										
Manual Roll Cutter & Holder 50cm A070 GM632CM	Technical Data					Sı	uitak	ole N	1ode	ls	
	This double-deck Manual Roll Cutting Device ensures that rolls are always cut by hand in the same size and serially.	A105	A100	A095	A090	A085	A080				
	 Contents Can be easily mounted on a mobile table, table or A050 Stand, Rolls of different widths up to 40cm can be cut at the same time, It is double-decker The roll slot has plastic bearings for easy roll rotation. Its double-sided, unbreakable blade is long-lasting, Upon request (optional) Side Shelf Set with Wall Mounting Equipment (A075) can be mounted on both sides of this stand. Pen holder and Indicator drawer are included in the Side Shelf Set. Thanks to the wall hanging equipment, it can be mounted on the wall and frees up a lot of space in the work area. All made of AISI 304 stainless steel. Dimension (L x W x H): 513 x 300 x 323 mm Weight: 7 Kg 										

Wall Hanger and Rack-Cup Set for A070 A075 GM632CMA	Technical Data					Sı	uita	ble N	1od	els		
	This Side Shelves and Hanging System is manufactured for the A070 Double-deck Manual Roll Cutting Apparatus.	A070										
	• 2 Metal Shelves + 1 Metal Cup + 1 Plastic Drawer and Metal Wall Hanger. • Side Shelf Set with Wall Mounting Equipment A070 can be mounted on either side of the Manual Cutting Attachment. Pen holder and Indicator drawer are included in the Side Shelf Set. • Manual Cutting Apparatus can be mounted on the wall thanks to its wall hanging equipment and frees up a lot of space in the work area. • Weight: 1.24 Kg	Y										
Manual Roll Cutter 75cm A065 GM632CC	Technical Data					Sı	uita	ble N	1od	els		
	This single-layer Manual Roll Cutting Apparatus ensures that rolls are always cut by hand in the same size and serially.	A105	A100	A095	A090	A085	A080	A077	A055	A050		
COUNTEN CONTRACTOR OF THE PROPERTY OF THE PROP	 Contents Can be easily mounted on a mobile table, table or A055 Stand, Rolls of different widths up to 60cm can be cut at the same time, The roll slot has plastic bearings for easy roll rotation. Its double-sided, unbreakable blade is long-lasting, Upon request (optional) Side Shelf Set (A077) can be mounted on both sides of this stand. All made of AISI 304 stainless steel. Dimension (L x W x H): 758 x 300 x 119 mm Weight: 8 Kg 						-					
Manual Roll Cutter 50cm A060 GM631CC	Technical Data					Sı	uita	ble N	1od	els		
	This single-layer Manual Roll Cutting Apparatus ensures that rolls are always cut by hand in the same size and serially.	A105	A100	A095	A090	A085	A080	A077	A055	A050		
GÜVOZ M	 Contents Can be easily mounted on a mobile table, table or A050 Stand, Rolls of different widths up to 40cm can be cut at the same time, The roll slot has plastic bearings for easy roll rotation. Its double-sided, unbreakable blade is long-lasting, Upon request (optional) Side Shelf Set (A077) can be mounted on both sides of this stand. Pen holder and Indicator drawer are included in the Side Shelf Set. All made of AISI 304 stainless steel. Dimension (L x W x H): 508 x 300 x 119 mm Weight: 4.80 Kg 											

Stand Wide 75cm A055 GM632CCS	Technical Data					Su	iitak	ole N	1ode	els				
	Stands allow the Rotary Sealers with Printer and Automatic Impulse Pouch Makers or cutting devices to be stacked on top of each other to free up space on the work table.	W509	20GM	46GM	45GM	40GM	30GM	20GM	10GM	M960	08GM	+W590	A065	A060
	 Contents Provides a very large work area on the table Can be easily mounted on a mobile table or desk, 30GM, 40GM Pouch Makers and A065 model Manual Roll Cutting Device can be mounted on this stand. 10GM and 20GM model Rotary Continuous Sealers with printer can be placed under the stand, the device can work in this case. Upon request (optional) Side Shelf Set (A077) can be mounted on both sides of this stand. Pen holder and Indicator drawer are included in the Side Shelf Set. The stand is stainless steel oven painted and it upper part is AISI 304 stainless steel. Dimension (L x W x H): 750 x 303 x 248 mm Weight: 5 Kg 													
Stand Narrow 50cm A050 GM631CCS	Technical Data					Su	iitak	ole N	1ode	els				
	Stands allow the Rotary Sealers without Printer and cutting devices to be stacked on top of each other to free up space on the work table.	W509	20GM	46GM	45GM	40GM	30GM	20GM	10GM	09GM	08GM	06GM+	05GD	A060
	 Contents Provides a very large work area on the table Can be easily mounted on a mobile table or desk, A060 model Manual Roll Cutting Device can be mounted on this stand. O5GD, 08GM and O9GM model Rotary Continuous Sealers with printer can be placed under the stand, the device can work in this case. Upon request (optional) Side Shelf Set (A077) can be mounted on both sides of this stand. Pen holder and Indicator drawer are included in the Side Shelf Set. The stand is stainless steel oven painted and it upper part is AISI 304 stainless steel. Dimension (L x W x H): 500 x 303 x 248 mm Weight: 3.65 Kg 													

Roller Conveyor XL 111cm A045 GM140RTXL	Technical Data					Su	ıitak	ole N	4od	els				
	This roller conveyor is ideal for preventing the packages from being damaged or torn by rubbing on the table while packaging the products in the devices. At the same time, there is no need for personnel to hold the packages along the device as the packages move through the device.	ВМ	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W960	W580	+W590	05GD	04GD+
	 Contents It has a stainless steel case and easy-rolling plastic bearings. Dimension (Length x Width x Height): 1110 x 260 x 50 / 150mm (Height Adjustable) Weight: 6.68 Kg Manufactured for 45GM, 46GM, 50GM and 60GM model Automatic Rotary Pouch Makers. 													
Roller Conveyor L 73cm A040 GM140RTL	Technical Data					Su	ıitak	ole N	Mode	els				
	This roller conveyor is ideal for preventing the packages from being damaged or torn by rubbing on the table while packaging the products in the devices. At the same time, there is no need for personnel to hold the packages along the device as the packages move through the device.	MDC	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W560	W580	+W590	05GD	04GD+
	 Contents It has a stainless steel case and easy-rolling plastic bearings. Dimension (Length x Width x Height): 760 x 260 x 50 / 150mm (Height Adjustable) Weight: 4.54 Kg Manufactured for 10GM and 20GM model devices 	_			_									
Roller Conveyor S 59cm A035 GM140RT	Technical Data					Su	ıitak	ole N	4od	els				
	This roller conveyor is ideal for preventing the packages from being damaged or torn by rubbing on the table while packaging the products in the devices. At the same time, there is no need for personnel to hold the packages along the device as the packages move through the device.	0GM	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W960	W580	+W590	05GD	04GD+
	 Contents It has a stainless steel case and easy-rolling plastic bearings. Dimension (Length x Width x Height): 560 x 260 x 50 / 150mm (Height Adjustable) Weight: 3.60 Kg Manufactured for 05GD, 08GM and 09GM model devices. 													

Side Rack Sets A030 GM14WTA for Rotary Sealers	Technical Data					Su	ıitak	ole N	1ode	els				
	These side racks are designed for Rotary Continuous Sealers . Thanks to the drawer and pen holder they contain, you can use indicators, pens, etc. used at CSSD. Ideal for storing ingredients	W509	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W560	W580	+W590	05GD	04GD+
	 Contents Can be mounted on both sides of Automatic rotary continuous sealers. Pen Holder and Indicator Drawer are included in the 3 pieces of metal Side Shelf It is AISI 304 stainless steel. Dimension (Length x Width x Height): 250 x 200 x 150 mm Weight: 1.00 Kg Manufactured for 10GM, 20GM, 45GM, 46GM, 50GM, 60GM model sealers. 													
Side Rack Sets A077 GM631CCA for Roll Cutter & Stands	Technical Data					Su	ıitak	ole N	1ode	els				
	These side racks are designed for Manual Roll Cutter & Stands. Thanks to the drawer and pen holder they contain, you can use indicators, pens, etc. used at CSSD. Ideal for storing ingredients	A070	A066	A065	A055	A050								
	 Contents Can be mounted on both sides of Manual Roll Cutter & Stands . Pen Holder and Indicator Drawer are included in the 4 pieces of metal Side Shelf It is AISI 304 stainless steel. Dimension (Length x Width x Height): 250 x 200 x 150 mm Weight: 1.20 Kg Manufactured for A050 and A055 model Stands, and A060, A065 and A070 model Manual Roll Cutters 													
Sorting Table - for Automatic Impulse Pouch Makers A025 GM450TC	Technical Data					Su	ıitak	ole N	1ode	els				
	This grouping tray is ideal for collecting the produced pouches in the same place and not scattering them when thrown from the device.	W509	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W560	W580	+W590	05GD	04GD+
Tour man and the state of the s	 Contents The cut pouches do not mix together. It is AISI 304 stainless steel. Dimension (Length x Width x Height): 563 x 274 x 67 mm Weight: 3.48 Kg Manufactured for 28GM, 30GM and 40GM model Automatic Impulse Pouch Makers 													

Eront Toble VI Go Dotom Law ith winter														
Front Table XL - for Rotary sealers with printer A020 GM111TC	Technical Data					Su	ıitak	ole N	1ode	els				
	The front table helps the user in the packaging process and prevents package damage.	W509	50GM	46GM	45GM	40GM	30GM	20GM	10GM	W960	W580	+W590	05GD	04GD+
Cannos and the second s	 Contents It is AISI 304 stainless steel. Dimension (Length x Width x Height): 540 x 305 x 81 mm Weight: 2.52 Kg Manufactured for 10GM and 20GM model Rotary Continuous Sealers with printer. 													
Front Table L for 06GM+ sealers A018 GM040TC	Technical Data					Su	ıitak	ole N	1od	els				
	The front table helps the user in the packaging process and prevents package damage.	W509	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W560	08GM	+W590	05GD	04GD+
	 Contents It is AISI 304 stainless steel. Dimension (Length x Width x Height): 535 x 318 x 69 mm Weight: 2.30 Kg Manufactured for O6GM+ Manual Impulse Pouch Maker 													
Front Table S - for narrow sealers A010 GD030TC	Technical Data					Su	ıitak	ole N	1od	els				
	The front table helps the user in the packaging process and prevents package damage.	W509	50GM	46GM	45GM	40GM	30GM	20GM	10GM	W960	08GM	+W590	05GD	04GD+
Concent.	 Contents It is AISI 304 stainless steel. Dimension (Length x Width x Height): 418 x 203 x 69 mm Weight: 1.58 Kg Manufactured for 02GD+, 03GD+, 04GD+, 05GD, 08GM and 09GM Sealers. 													

Reel holder Stand Large A016 GM040RC	Technical Data	Suitable Models																			
	When the device is hung on the wall to benefit from the working area, a Wall Hanging Apparatus is needed so that the rolls can stand on the device.	W509	20GM	46GM	45GM	40GM	30GM	20GM	10GM	WD60	W580	+W590	05GD	04GD+							
	 Contents It is AISI 304 stainless steel. Dimension (Length x Width x Height): 535 x 180 x 72 mm Weight: 1.50 Kg Manufactured for O6GM+ Manual Impulse Pouch Maker 																				
Reel holder Stand Small A015 GD030RC	Technical Data	Suitable Models											Suitable Models								
	When the device is hung on the wall to benefit from the working area, a Wall Hanging Apparatus is needed so that the rolls can stand on the device.	W509	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W960	W580	+W590	05GD	04GD+							
	 Contents It is AISI 304 stainless steel. Dimension (Length x Width x Height): 418 x 180 x 72 mm Weight: 1.00 Kg Manufactured for 02GD+, 03GD+ and 04GD+ model Manual Impulse Pouch Makers. 																				
ABGPack® Medical Packaging - Sterilization Reels A999 Medipack	Technical Data	Suitable Models											Suitable Models								
Rulo MEDIKAL STERILIZASYON RULOSU Q	ABGPack® sterilization rolls are produced under EN ISO 9001:2015 & ISO 13485 quality systems and comply with EN868-5	W509	20GM	46GM	45GM	40GM	30GM	20GM	10GM	W560	W580	+W590	05GD	04GD+							
	 Contents Solvent-free ABGPack® sterilization rolls have excellent adhesion ability, ability to be stable against sterilization effects, excellent mechanical strength and excellent barrier properties in the same strong multi-layered structure. With its solvent-free structure, it provides full protection from solvent migration and contamination risks. With ABGPack®, the weld quality in adhesion areas can be easily determined and offers the advantage of being used safely at every stage of sterilization processes. 																				


ABGIndicators® Sterilization Monitoring Indicators A997 GM997STM	Technical Data	Suitable Models										
	ABGIndicators® & Biological and Chemical Indicators & Labels and Documentation System with/without indicators. For Steam, Ethylene Oxide, Formaldehyde and Hydrogen Peroxide / Plasma sterilization processes	Steam	EO	FORM	H202							
	 Contents The sterilization documentation system is used for patient-related traceability of all medical sterilization processes. The documentation page contains all the information necessary to release the sterilization batch. One sterilizer and one documentation sheet for the day are used. Using a labeling device (Label Gun), the documentation sheet and each package can be labeled with a self-adhesive label containing the following information: Production and expiration date, Person in charge, Contents of the packaging, Sterilizer and batch number. 											
ABGTape® Identification Tapes for Surgical Instruments A996 GM996IDT	Technical Data	Suitable Models										
	It provides organization and tracking of your valuable tools easily and economically. It is resistant to Steam, Ethylene Oxide, Formaldehyde and Hydrogen Peroxide / Plasma sterilization processes.	Steam	EO	FORM	H202							
	 Contents Surgical Instrument Identification Tapes adhere securely, do not break or flake, and do not come off or unravel on their own. Available in 275 colors and styles (symbols and characters), color-coding options are nearly endless. When replacement is required, it can be easily removed and a new color and style combination can be made with the tool and "Adhesive Remover", which is environmentally friendly, does not contain toxic solvents, has no harmful smell and is compatible with many surfaces. The products do not contain Latex and are not sterile The product is produced in the form of 312, 352, 376 and 400 cut strips on an A4 page, in two main groups: 1- Plain and/or Patterned Tapes and 2- Lettered and/or Numbered. There are options with different colors and symbols for all of these groups. 											

Seal-Check-Test - Tyvek® A170 GM170SCT	Technical Data	Suitable Models												
	This test is produced to test whether the pouches produced for H2O2/Plasma Sterilization are working in the sealer, at factory settings and in accordance with the standard defined in EN-ISO 11607-2, before the first use every day.		20GM	46GM	45GM	40GM	30GM	20GM	10GM	W560	W580	+W590	05GD	04GD+
SEAL CHECK-TYVEK® Testing for sealing machines Kapatma Chazlar testi igin Please put this lates sheet ire to a pouch and pass through the sealer cnce. The seal must be darker and homogeneous than the environment throughout. Sealer / Chaz: Location / Yer: Approval / Onay: Date / Tarih: REF: 222-110-070 Tel: +90 212 13 5353 per bar. +90 212 212 2498 Luflen bu lest kajd-din popelinize yerinsz. per kar. +90 212 212 2498 Yap-shrmak karminn diger kusmidan daha koyu ve homogen olmasa gerekmektedir. Sealer / Chaz: Date / Tarih: Remarks / Notiar:	 Contents At the bottom of the Test Sheet, there are sections containing the serial number of the device, the place where the device is used, the person performing the test and his/her signature, and comment/approval areas about the test result. One package contains 100 test sheets. Dimensions 175 x 75 x 15 mm [Length x Width x Height] Weight 0.180Kg. 													
Seal-Check-Test - Standard A165 GM165SCS	Technical Data		Suitable Models											
	This test is produced to test whether the pouches produced for STEAM / EO / FORM Sterilization are working in the sealer, at factory settings and in accordance with the standard defined in EN-ISO 11607-2, before the first use every day.	WD09	20GM	46GM	45GM	40GM	30GM	20GM	10GM	WD60	08GM	+W590	05GD	04GD+
SEAL CHECK Testing for sealing machines Kapatma Cihazları testi için Please put this test sheet into a pouch and pass through the sealer once. The seal must be darker and homogeneous than the environment throughout. Luften bu test Kagdum poşehitüre yerleşting haza veriniz. Yapıştırma kısımının diğer kısımdan daha koyu ve homojen olması gerekmektedir. Sealer / Cihaz: REF: 222-110-060 LOT: 1000160515 REF: 222-1110-060 Cikiden Food Pood Pood Pood Pood Pood Pood Pood	 Contents At the bottom of the Test Sheet, there are sections containing the serial number of the device, the place where the device is used, the person performing the test and his/her signature, and comment/approval areas about the test result. One package contains 250 test sheets. Dimensions 175 x 75 x 22 mm [Length x Width x Height] Weight 0.195 Kg 													

MANUFACTURER COMPANY SPECIALIZED IN MEDICAL SEALERS SINCE 1990.

Our medical packaging devices and accessories meet the requirements of the World Health Organization (WHO) guideline as well as the current guideline of the World Federation of Hospital Sterilization Sciences (WFHSS).

It ensures validable sealing according to EN ISO 11607-2 and the international guideline ISO/TS 16775. Document this in writing with a Validation Report every time it is used for the first time or whenever requested.



REGISTERED TRADEMARKS OF GÜNDEM A.Ş.

GÜNDEM® | ABGMat® | ABGPack® | ABGTape® | ABGIndicators®

FOR OTHER PRODUCTS AND MODELS: www.gundemas.com





Gündem produces the following products, which are unique and under international patent decision, in its own facilities, developed with the experience it has gained since 1990. Our medical packaging devices and accessories meet the requirements of the World Health Organization (WHO) guideline as well as the current guideline of the World Federation of Hospital Sterilization Sciences (WFHSS). It ensures validable sealing according to EN ISO 11607-2 and the international guideline ISO/TS 16775.

Document this in writing with a Validation Report every time it is used for the first time or whenever requested.

Our other productions

* Surgical Instrument Identification Tapes

* Sterilization Packaging Products

* Chemical and Biological Indicators and accessories for sterilization monitoring

Address: Mecidiyekoy Mah., Sehit Er Cihan Namli Cad., Mert Sok., No: 1/A 34381 Sisli - Istanbul - Türkiye Tel.: +90 212 213 5353 Pbx., Fax.: +90 212 212 2498 E-Mail: order@gundemas.com Web.: www.gundemas.com



Quality

Management

System